

2018:

Albanna W, Kotliar K., Lüke JN, Alpdogan S, Conzen C, Lindauer U, Clusmann H, Hescheler J, Vilser W, **Schneider T**, Schubert GA (2018) Non-invasive evaluation of neurovascular coupling in the murine retina by dynamic retinal vessel analysis. PLoS One 13 (10) doi: 10.1371/journal.pone.0204689. eCollection 2018.

Neumaier Felix, Alpdogan Serdar, Hescheler Jürgen, **Schneider Toni** (2018) A practical guide to the preparation and use of metal ion-buffered systems in physiological research. Acta Physiologica (Oxf.) 222(3). doi: 10.1111/apha.12988. [Epub ahead of print] Review.

Neumaier Felix, Alpdogan Serdar, Hescheler Jürgen, **Schneider Toni** (2018) Protein phosphorylation maintains the normal function of cloned human Ca_v2.3 channels. J. Gen. Physiology, in press. doi: 10.1085/jgp.201711880.

Neumaier Felix, Paterno Mario, Alpdogan Serdar, Tevoufouet Etienne E., **Schneider Toni**, Hescheler Jürgen, Albanna Walid (2018) In Reply to „Corpus Callosotomy for Drug-Resistant Schizophrenia: Novel Treatment Based on Pathophysiology”. World Neurosurgery, in press, doi: 10.1016/j.wneu.2018.04.206.

Walid Albanna*; Felix Neumaier*; Jan Niklas Lüke; Konstantin Kotliar; Catharina Conzen; Ute Lindauer; Jürgen Hescheler; Hans Clusmann; **Toni Schneider**[§]; Gerrit Alexander Schubert[§] (2018) Unconjugated bilirubin modulates neuronal signaling only in wildtype mice, but not after ablation of the R-type / Ca_v2.3 voltage-gated calcium channel. CNS Neuroscience & Therapeutics, 24(3):222-230. doi: 10.1111/cns.12791. */[§] Contributed equally

Albanna W., Lüke J.N., Neumaier F., Alpdogan S., Conzen, C., Hescheler J., Clusmann H., Lüke M., Schubert G.A., **Schneider T**. (2018) Selected aspects of retinal signaling and energy metabolism and its perspective as a cerebral surrogate model. New Frontiers in Ophthalmology, in press.

Neumaier F., Akhtar I, Hescheler J., **Schneider T**. (2018) Reciprocal modulation of Ca_v2.3 E-/ R-Type voltage-gated calcium channels by copper(II) ions and kainate. J. Neurochem., in press, doi: 10.1111/jnc.14546.

Schneider T, Alpdogan S., Hescheler J., Neumaier F. (2018) In vitro and in vivo phosphorylation of the Ca_v2.3 voltage-gated R-type calcium channel. Channels, in press.

2017:

Neumaier Felix, Mario Paterno, Serdar Alpdogan, Etienne E. Tevoufouet, **Toni Schneider**, Jürgen Hescheler, Walid Albanna (2017) Surgical approaches in psychiatry: A survey of the world literature on psychosurgery. World Neurosurgery, 97: 140-147. doi: 10.1016/j.wneu.2016.10.008. Review.

Tevoufouet E.E., Nembo E.N., Distler F., Neumaier F., Hescheler J., Nguemo F., **Schneider T**. (2017) Multiple Nickel-Sensitive Targets Elicit Cardiac Arrhythmia in Isolated Mouse Hearts after Pituitary Adenylate Cyclase-Activating Polypeptide-Mediated Chronotropy. Pharmacological Research 117, 140-147.

Marcel A. Kamp, Jasper v. Lieshout, Maxine Dibué-Adjei, Jasmin Weber, **Toni Schneider**, Tanja Restin, Igor Fischer, Hans-Jakob Steiger (2016) A systematic and meta-analysis of mortality in experimental mouse models analyzing delayed cerebral ischemia after subarachnoid hemorrhage. Translational Stroke Research, 8: 206 – 219. doi: 10.1007/s12975-016-0513-3.

Mahdy Ranjbar, **Toni Schneider**, Salvatore Grisnati, Carl Brand, Julia Lüke, Matthias Lüke (2017) The effect of anakinra on retinal function in isolated perfused vertebrate retina. J. Curr. Ophthalmology, 29:69-71. doi: 10.1016/j.joco.2016.12.002. eCollection 2017 Mar.

Poomvanicha Montatip, Matthes Jan, Domes Katrin, Patrucco Enrico, Angermeier Elisabeth, Laugwitz Karl-Ludwig, **Schneider Toni**, Hofmann Franz (2017) Beta-adrenergic regulation of the heart expressing the Ser1700A/Thr1704A mutated Ca_v1.2 channel. Molecular and Cellular Cardiology 111: 10-16.

Christian Ritter*, Martin K. R. Svačina*, Ilja Bobylev, Abhijeet Joshi, **Toni Schneider**, Helmar C. Lehmann (2017) Impact of age and polytherapy on fingolimod induced bradycardia: a preclinical study. *Journal of Neuroimmune Pharmacology*, 12: 204-209. * contributed equally

Jasper van Lieshout, Maxine Dibué-Adjei, Jan F. Cornelius, Philipp J. Slotty, **Toni Schneider**, Restin T., Boongaaarts H.D., Hans-Jakob Steiger, Athanasios K. Petridis, Marcel A. Kamp (2017) An Introduction to the pathophysiology of aneurysmal subarachnoid hemorrhage. *Neurosurgical Review*, doi: 10.1007/s10143-017-0827-y. [Epub ahead of print].

Walid Albanna, Miriam Weiss, Catharina Conzen, Hans Clusmann, **Schneider Toni**, Nicole Heussen, Martin Reinsch, Marguerite Müller, Anke Höllig, Gerrit Alexander Schubert (2017) Systemic and cerebral concentration of nimodipine during established and experimental vasospasm treatment. *World Neurosurgery*, 102: 459-465.

Walid Albanna, Jan Niklas Lueke, Volha Sjapic, Konstantin Kotlair, Jürgen Hescheler, Hans Clusmann, Sergej Sjapic, **Toni Schneider**, Gerrit Alexander Schubert, Felix Neumaier (2017) Electroretinographic assessment of inner retinal signaling in the isolated and superfused murine retina. *Current Eye Research*, 42: 1518-1526.

Dibué-Adjei Maxine, Marcel A. Kamp, Serdar Alpdogan, Etienne E. Tevoufouet, Wolfram F. Neiss, Jürgen Hescheler, **Toni Schneider** (2017) Ca_v2.3 (R-Type) Calcium Channels are Critical for Mediating Anticonvulsive and Neuroprotective Properties of Lamotrigine In Vivo. *Cell. Physiol. Biochem.*, 44: 935 – 947.

Schneider T. (2017) Response to “Radiotelemetric EEG recordings in mice — Critical role of transmitter bandwidth and nominal sampling rate in frequency analysis”. *Epilepsy Res.*, 132: 117-118. <http://dx.doi.org/10.1016/j.eplepsyres.2016.12.004>

2015:

Schneider T, Dibue-Adjei M (2015) Ca_v2.3 E-R-Type Voltage-Gated Calcium Channels Modulate Sleep in Mice. *Sleep* 38:499.

Drobinskaya I, Neumaier F, Pereverzev A, Hescheler J, **Schneider T** (2015) Diethyldithiocarbamate-mediated zinc ion chelation reveals role of Cav2.3 channels in glucagon secretion. *Biochim Biophys Acta Mol Cell Res* 1854:953-964.

Ranjbar M, **Schneider T**, Brand C, Grisanti S, Lüke J, Lüke M (2015) The effect of disease-modifying anti-rheumatic drugs on retinal function in the electrophysiological ex-vivo model of the isolated perfused vertebrate retina. *Ophthalmic Res* 53:136-140.

Neumaier F, Dibué-Adjei M, Hescheler J, **Schneider T** (2015) Voltage-gated calcium channels: Determinants of channel function and modulation by inorganic cations. *Progress in Neurobiology* 129:1-36.

Kamp MA, Dibué M, Sommer C, Steiger HJ, **Schneider T**, Hänggi D (2015) Evaluation of a murine single-blood-injection SAH model. *PLoS One* 9:e114946.

Alnawaiseh M, Albanna W, Abumuaileq R, Böhm MRR, Eter N, **Schneider T** (2015) Einfluss von Hypericin auf die isolierte pigmentepithelfreie Netzhaut: Eine elektroretinographische Untersuchung. *Klinische Monatsblätter für Augenheilkunde* 232: 1304-1307..

Siapich S, Akhtar I, Hescheler J, **Schneider T**, Lüke M (2015) Low concentrations of ethanol but not of dimethyl sulfoxide (DMSO) impair reciprocal retinal signal transduction. *Graefes Arch Clin Exp Ophthalmol* 253:1713-1719.

Schneider T, Dibué-Adjei M, Neumaier F, Akhtar I, Hescheler J, Kamp MA, Tevoufouet EE (2015) R-type voltage-gated Ca²⁺ channels in cardiac and neuronal rhythmogenesis. *Current Molecular Pharmacology* 8:102-108.

2014:

Kamp MA, Dibué M, Sommer C, Steiger HJ, **Schneider T**, Hänggi D (2014) Evaluation of a murine single-blood-injection SAH model. *PLoS One* 29: e114946. Doi: 10.1371/journal.pone.0114946. eCollection 2014..

Tevoufouet EE, Nembo EN, Dibue-Adjei M, Hescheler J, Nguemo F, **Schneider T** (2014) Cardiac Functions of Voltage-Gated Ca(2+) Channels: Role of the Pharmacoresistant Type (E-/R-Type) in Cardiac Modulation and Putative Implication in Sudden Unexpected Death in Epilepsy (SUDEP). *Rev Physiol Biochem Pharmacol* 167: 115-139.

Schneider T, Dibue-Adjei M (2014) Ca_v2.3 E-/R-Type Voltage-Gated Calcium Channels Modulate Sleep in Mice. *Sleep* Epub ahead of print (Letter to the editor).

Dibué M, Kamp MA, Neumaier F, Steiger H-J, Hänggi D, Hescheler J, **Schneider T** (2014) Cardiac Phenomena During Kainic-acid Induced Epilepsy and Lamotrigine Antiepileptic Therapy. *Epilepsy Res* 108: 666-674.

Tura A, Alt A, Haritoglou C, Meyer CH, **Schneider T**, Grisanti S, Lüke J, Lüke M, for the international Chromovitrectomy Collaboration (2014) Testing the effects of the dye Acid violet-17 on retinal function for an intraocular application in vitreo-retinal surgery. *Graefes Arch Clin Exp Ophthalmol* 252: 1927-1937.

Ranjbar M, Alt A, Nassar K, Reinsberg M, **Schneider T**, Grisanti S, Lüke J, Lüke M (2014) The Concentration-Dependent Effects of Indocyanine Green on Retinal Function in the Electrophysiological ex vivo Model of Isolated Perfused Vertebrate Retina. *Ophthalmic Res* 51: 167-171.

2013:

Kamp MA, Dibué MAH, Etminan N, Steiger H-J, **Schneider T**, Hänggi D (2013) Evidence for direct impairment of neuronal function by subarachnoid metabolites following SAH. *Acta Neurochirurgica* 155: 255-260

Alt A, Hilgers R-D, Tura A, Nassar K, **Schneider T**, Hüber A, Januschowski K, Grisanti S, Lüke J, Lüke M (2013) The Neuroprotective Potential of Rho-Kinase Inhibition in Promoting Cell Survival and Reducing Reactive Gliosis in Response to Hypoxia in Isolated Bovine Retina. *Cell Physiol Biochem* 32: 218 - 234.

Dartsch T, Fischer R, Gapelyuk A, Weiergräber M, Ladage D, **Schneider T**, Schirdewan A, Reuter H, Mueller-Ehmsen J, Zobel C (2013) Aldosterone induces electrical remodeling independent of hypertension. *Int J Cardiol* 164: 170-178.

Dibué M, Tevoufouet EE, Neumaier F, Krieger A, Kiel A, Evdokimov D, Galetin T, Alpdogan S, Akthar I, Scharf S, Clemens R, Radhakrishnan K, Hescheler J, **Schneider T**, Kamp MA (2013) Protein interaction partners of Cav2.3 R-type voltage-gated calcium channels. In: *Presynaptic Voltage-gated Calcium Channels* (Stephens GJ, Mochida S, eds; Basel, Switzerland: MDPI), Pharmaceuticals 6: 151-174.

Dibué MAH, Alpdogan S, Tevoufouet EE, Hescheler J, Kamp M, **Schneider T** (2013) Cav2.3 (R-type) Calcium Channels are critical for Mediating Anticonvulsive and Neuroprotective Properties of Lamotrigine and Topiramate in vivo. *Epilepsia* 54: 1542-1550.

Galetin T, Tevoufouet EE, Sandmeyer J, Matthes J, Nguemo F, Hescheler J, Weiergräber M, **Schneider T** (2013) Pharmacoresistant Ca_v2.3 (E-/R-type) voltage-gated calcium channels influence heart rate dynamics and contribute to cardiac impulse conduction. *Cell Biochemistry & Function* 31: 434 – 449.

Schneider T, Dibué M, Hescheler J (2013) How "Pharmacoresistant" is Cav2.3, the Major Component of Voltage-gated R-type Ca²⁺ Channels? In: *Calcium Antagonists* (**Schneider T**, Stephens GJ, eds; Basel, Switzerland: MDPI), Pharmaceuticals 6: 759-776.

Rassaei M, Thelen M, Abumuaileq R, Hescheler J, Luke M, **Schneider T** (2013) Effect of high-intensity irradiation from dental photopolymerization on the isolated and superfused vertebrate retina. *Graefes Arch Clin Exp Ophthalmol* 251: 751-762.

Yang L, Topia I, **Schneider T**, Stephens GJ (2013) Phorbol ester modulation of Ca²⁺ channels mediate nociceptive transmission in dorsal horn neurones. In: Calcium Antagonists (**Schneider T**, Stephens GJ, eds; Basel, Switzerland: MDPI), Pharmaceuticals 6: 777-787.

2012:

Kamp MA, Hanggi D, Steiger HJ, **Schneider T** (2012) Diversity of presynaptic calcium channels displaying different synaptic properties. *Rev Neurosci* 23: 179-190.

Kamp MA, Shakeri B, Tevoufouet EE, Krieger A, Henry M, Behnke K, Herzig S, Hescheler J, Radhakrishnan K, Parent L, **Schneider T** (2012) The C-terminus of human Ca(v)2.3 voltage-gated calcium channel interacts with alternatively spliced calmodulin-2 expressed in two human cell lines. *Biochim Biophys Acta* 1824: 1045-1057.

Müller R, Struck H, Ho MS, Brockhaus-Dumke A, Klosterkötter J, Broich K, Hescheler J, **Schneider T**, Weiergräber M (2012) Atropine-sensitive hippocampal theta oscillations are mediated by Ca(v)2.3 R-type Ca(2)(+) channels. *Neuroscience* 205: 125-139.

Lammers C, Dartsch T, Brandt MC, Rottlander D, Halbach M, Peinkofer G, Ockenpoehler S, Weiergräber M, **Schneider T**, Reuter H, Müller-Ehmsen J, Hescheler J, Hoppe UC, Zobel C (2012) Spironolactone prevents aldosterone induced increased duration of atrial fibrillation in rat. *Cell Physiol Biochem* 29: 833-840.

Parajuli LK, Nakajima C, Kulik A, Matsui K, **Schneider T**, Shigemoto R, Fukazawa Y (2012) Quantitative Regional and Ultrastructural Localization of the Cav2.3 Subunit of R-type Calcium Channels in Mouse Brain. *J Neurosci* 32: 13555-13567.

Kamp MA, Dibué MAH, **Schneider T**, Steiger H-J, Hänggi D (2012) Calcium and Potassium Channels in Experimental Subarachnoid Hemorrhage and Transient Global Ischemia. *Stroke Research and Treatment* doi: 10.1155/2012/382146.

Kamp MA, Dibué MAH, Etminan N, Steiger H-J, **Schneider T**, Hänggi D (2012) Evidence for direct impairment of neuronal function by subarachnoid metabolites following SAH. *Acta Neurochirurgica DOI* 10.1007/s00701-012-1559-y.

2011:

Alnawaiseh M, Albanna W, Chen C-C, Campbell KP, Hescheler J, Lüke M, **Schneider T** (2011) Two separate Ni²⁺ sensitive voltage-gated Ca²⁺ channels modulate transretinal signalling in the isolated murine retina. *Acta Ophthalmologica* 89: e579-e590.

Radhakrishnan K, Krieger A, Maxine Dibué, Hescheler J, **Schneider T** (2011) APLP1 and Rab5A Interact with the II-III loop of the Voltage-gated Ca²⁺-channel Cav2.3 and Modulate its Internalization Differently. *Cellul Physiol Biochem* 28: 603-612.

Radhakrishnan K, Kamp MA, Siapich SA, Hescheler J, Luke M, **Schneider T** (2011) Ca(v)2.3 Ca²⁺ Channel Interacts with the G1-subunit of V-ATPase. *Cell Physiol Biochem* 27: 421-432.

Radhakrishnan K, Krieger A, Tevoufouet EE, Bähr M, Hescheler J, **Schneider T** (2011) Interaction of Recombinant and Native Cav2.3 E-/R-type Voltage-gated Ca²⁺ Channels with the Molecular Chaperone Hsp70. *J Receptor Ligand Channel Research* 4: 29-40.

Alnawaiseh M, Albanna W, Banat M, Abumuailleq R, Hescheler J, **Schneider T** (2011) Electroretinographic Recordings from the Isolated and Superfused Murine Retina. In: *Electroretinograms* (Belulsic G, ed), pp 175-190. Rijeka: Croatia.

2010:

Bourdin B, Marger F, Wall-Lacelle S, **Schneider T**, Klein H, Sauve R, Parent L (2010) Molecular determinants of the CaV β -induced plasma membrane targeting of the CaV1.2 channel. *J Biol Chem* 285: 22853-22863.

Galetin T., Weiergräber M., Hescheler J., **Schneider T.** (2010) Analyzing murine electrocardiogram with PhysioToolkit. *J. Electrocardiol.*, 43: 701-705.

Lee S, Grafweg S, **Schneider T**, Jimenez M, Giacobino JP, Ghanem A, Tiemann K, Bloch W, Müller-Ehmsen J, Schwinger RH, Brixius K (2010) Total beta-adrenoceptor deficiency results in cardiac hypotrophy and negative inotropy. *Physiol Res* 59: 679-689.

Lüke M, Januschowski K, Tura A, Lüke J, Nassar K, Lüke C, **Schneider T**, Szurman P, Grisanti S, Bartz-Schmidt KU (2010) Effects of pegaptanib sodium on retinal function in isolated perfused vertebrate retina. *Curr Eye Res.* 35:248-254.

Radhakrishnan K., Hescheler J., **Schneider T.** (2010) Heat shock proteins and ion channels. Functional interactions and therapeutic consequences. *Curr Pharm. Biotechnol.* 11: 175-179.

Siapich S.A., Wrubel H., Albanna W., Alnawaiseh M., Hescheler J., Weiergräber M., Lüke M., **Schneider T.** (2010) Effect of ZnCl₂ and chelation of zinc ions by N,N-diethyldithiocarbamate (DEDTC) on the ERG b-wave amplitude from the isolated superfused vertebrate retina. *Curr. Eye Res.* 35: 322-334.

Vasudevan A, Ho MS, Weiergräber M, Nischt R, **Schneider T**, Lie A, Smyth N, Köhling R. (2010) Basement membrane protein nidogen-1 shapes hippocampal synaptic plasticity and excitability. *Hippocampus.* 2010 May;20(5):608-20

2009:

Albanna W, Banat M, Albanna N, Alnawaiseh M, Siapich SA, Igelmund P, Weiergräber M, Lüke M, **Schneider T.** (2009) Longer lasting electroretinographic recordings from the isolated and superfused murine retina. *Graefes Arch Clin Exp Ophthalmol.* 247:1339-1352 (Note the animation in the supplement!).

Bartels P, Behnke K, Michels G, Groner F, **Schneider T**, Henry M, Barrett PQ, Kang HW, Lee JH, Wiesen MH, Matthes J, Herzig S (2009) Structural and biophysical determinants of single Ca(V)_{3.1} and Ca(V)_{3.2} T-type calcium channel inhibition by N(2)O. *Cell Calcium* 46: 293-302.

Hildebrand ME, Isope P, Miyazaki T, Nakaya T, Garcia E, Feltz A, **Schneider T**, Hescheler J, Kano M, Sakimura K, Watanabe M, Dieudonné S, Snutch TP (2009) Functional coupling between mGluR1 and Cav3.1 T-type calcium channels contributes to parallel fiber-induced fast calcium signaling within Purkinje cell dendritic spines. *J Neurosci.* 29:9668-82

Joksovic P.M., Weiergräber M., Lee W.Y., Struck H., **Schneider T.**, Todorovic S.M. (2008) Isoflurane-sensitive presynaptic R-type calcium channels contribute to inhibitory synaptic transmission in the rat thalamus. *J.Neurosci.*, 29: 1434-1445..

Matar Nora, Jin Wei, Wrubel Heiko, Hescheler Jürgen, **Schneider Toni**, Weiergräber, Marco (2009) Zonisamide Block of Cloned Human T-Type Voltage-Gated Calcium Channels. *Epilepsy Research*, 83: 224-234.

Siapich S.A., Banat M., Albanna W., Hescheler J., Lüke M., **Schneider T.** (2009) Antagonists of ionotropic GABA receptors impair the NiCl₂ mediated stimulation of the ERG b-wave amplitude from the isolated superfused vertebrate retina. *Acta Ophthalmol.*, 87: 854-865.

2008:

Salvador-Recatalà, Vicenta, **Schneider, Toni** and Robert M. Greenberg (2008) Atypical properties of a conventional Ca₂₊ channel β subunit from the platyhelminth *Schistosoma mansoni*. *BMC Physiol.*, 8: 6.

Matthias Lüke, Kai Januschowski, Julia Beutel, Max Warga, Salvatore Grisanti, Swaantje Peters, **Toni Schneider**, Christoph Lüke, Karl U Bartz-Schmidt, Peter Szurman (2008) The effects of triamcinolone crystals on retinal function in a model of isolated perfused vertebrate retina. *Experimental Eye Res.*, 87, 22-29.

Weiergräber, Marco, Margit Henry, Matthew S.P. Ho, Henrik Struck, Jürgen Hescheler, **Toni Schneider** (2008) Altered Thalamocortical Rhythmicity in Ca_v2.3-Deficient Mice. *Mol.Cellul.Neurosci.*, 39: 605 - 618.

Weiergräber M., Hescheler J., **Schneider T.** (2008) [Human calcium channelopathies: Voltage-gated Ca₂₊ channels in etiology, pathogenesis, and pharmacotherapy of neurologic disorders.] *Nervenarzt* 79: 426-36; German.

Banat M., Lüke M., Siapich S.A., Hescheler J., Weiergräber M., **Schneider T.** (2008) The dihydropyridine isradipine inhibits the murine but not the bovine a-wave response of the electroretinogram. *Acta Ophthalmol.*, 86: 676-82.

2007:

Lüke M., Szurman P., **Schneider T.**, Lüke C. (2007) The effects of the phosphodiesterase type V inhibitor sildenafil on human and bovine retinal function *in vitro*. *Graefes Archive for Clinical and Experimental Ophthalmology*, 245: 1211-1215.

Weiergräber M., Henry M., Radhakrishnan K., Hescheler J., **Schneider T.** (2007) Hippocampal seizure resistance and reduced neuronal excitotoxicity in mice lacking the Ca_v2.3 E/R-type voltage-gated calcium channel. *J. Neurophysiology*, 97, 3660 - 3669.

Zhang Q., Bengtsson M., Partridge C., Salehi A., Braun M., Cox R., Eliasson L., Johnson P.R., Renström E., **Schneider T.**, Berggren P.-O., Göpel S., Ashcroft F.M., Rorsman P. (2007) R-type Ca₂₊-channel-evoked CICR regulates glucose-induced somatostatin secretion. *Nature Cell Biol.*, 9, 453-60.

Matthias Lüke, Krott Ralf, Warga Max, Szurman Peter, Grisanti Salvatore, Bartz-Schmidt Karl Ulrich, **Schneider Toni**, Lüke Christoph (2007) Effects of the Protein Tyrosine Kinase Inhibitor Genistein and Taurine on Retinal Function in Isolated Superfused Retina. *Graefes Archive for Clinical and Experimental Ophthalmology*, 245, 242 – 248.

Lüke M, Januschowski K, Warga M, Beutel J, Leitritz M, Gelisken F, Grisanti S, **Schneider T**, Lüke C, Bartz-Schmidt KU, Szurman P (2007) The retinal tolerance to bevacizumab in co-application with a recombinant tissue plasminogen activator. *Br J Ophthalmol* 91: 1077-1082.

2006:

Andreas Krieger*, Kayalvizhi Radhakrishnan*, Alexey Pereverzev*, Siarhei A. Siapich, Mohammed Banat, Marcel A. Kamp, Jérôme Leroy, Udo Klöckner, Jürgen Hescheler, Marco Weiergräber, **Toni Schneider** (2006) The molecular chaperone hsp70 interacts with the cytosolic II-III loop of the Cav2.3 E-type voltage-gated Ca₂₊ channel. *Cell Physiol Biochem.* 17, 97-110; (* = contributed equally)

Lüke C., Lüke M., Sickel W., **Schneider T.** (2006) Effects of patent blue on human retinal function. *Graefes Archive for Clinical and Experimental Ophthalmology*, 244, 1188 - 90.

Köhling R., Nischt R., Vasudevan A., Ho M., Weiergräber M., **Schneider T.**, Smyth N. (2006) Neurodegener. Dis. 3, 56 – 61.

Weiergräber M., Kamp M.A., Radhakrishnan K., Hescheler J., **Schneider T.** (2006) The Cav2.3 voltage-gated Ca₂₊ channel in epileptogenesis – Shedding new light on an enigmatic channel. *Neurosci. Biobehav. Rev.*, 30, 1122-1144.

Natrajan R., Little S.E., Reis-Filho J.S., Hing L., Messahel B., Grundy P.E., Dome J.S., **Schneider T.**, Vujanic G.M., Pritchard-Jones K., Jones C. (2006) Amplification and overexpression of CACNA1E correlates with relapse in favorable histology Wilms' tumors. *Clin. Cancer Res.*, 12, 7284 – 93.

Lüke M., Warga M., Ziemssen F., Gelisken F., Grisanti S., **Schneider T.**, Lüke C., Partsch M., Bartz-Schmidt K.U., Szurman P., for the Tuebingen Bevacizumab Study Group (2006) Effects of bevacizumab on retinal function in isolated vertebrate retina. *British Journal of Ophthalmology*, 90, 1178 - 1182

Marco Weiergräber, Margit Henry, Andreas Krieger, Marcel Kamp, Kayalvizhi Radhakrishnan, Jürgen Hescheler, **Toni Schneider** (2006) Altered seizure susceptibility in mice lacking the Cav2.3 E-type Ca₂₊ channel. *Epilepsia.* 47, 839-50.

2005:

Matthias Lüke , Marco Weiergräber, Carl Brand, Siarhei A. Siapich, Mohammed Banat, Jürgen Hescheler, Christoph Lüke, **Toni Schneider** (2005) The isolated perfused bovine retina - a sensitive tool for pharmacological research on retinal function .*Brain Res Brain Res Protoc.* 16, 27-36.

Matthias Lüke, Christoph Lüke, J. Hescheler, **T. Schneider**, W. Sickel (2005) Effects of Phosphodiesterase Type 5 Inhibitor Sildenafil on Retinal Function in Isolated Superfused Retina. *J. Ocul. Pharmacol. & Ther.*,21: 305-14

Matthias Lüke, Margit Henry, Thea Lingohr, Mehran Maghsoodan, Jürgen Hescheler, Marco Weiergräber, Werner Sickel, **Toni Schneider** (2005) A Ni²⁺-sensitive component of the ERG-b-wave from the isolated bovine retina is related to E-type voltage-gated Ca₂₊ channels. *Graefes Archive for Clinical and Experimental Ophthalmology*, 243: 933-941.

Marco Weiergräber, Margit Henry, Jürgen Hescheler, Neil Smyth & **Toni Schneider** (2005) Electrographic and deep intracerebral EEG recording in mice using a telemetry system. *Brain Res. & Brain Res. Protoc.*, 14: 154-64.

Kamp M.A., Krieger A., Henry M., Hescheler J., Weiergräber M., **Schneider T.** (2005) Presynaptic "Ca_v2.3 containing" E-type Ca₂₊ channels share dual roles during neurotransmitter release. *Eur J Neurosci.* 21:1617-25. Review

Pereverzev A, Salehi A, Mikhna M, Renström E, Hescheler J, Weiergräber M, Smyth N, **Schneider T.** (2005) The ablation of the Ca_v2.3/E-type voltage-gated Ca₂₊ channel causes a mild phenotype despite an altered glucose induced glucagon response in isolated islets of Langerhans. *Eur. J. Pharmacol.* 511, 65-72.

Xingjun Jing, Dai-Qing Li, Charlotta S. Olofsson, Albert Salehi, Vikas V. Surve, José Caballero, Rosita Ivarsson, Ingmar Lundquist, Alexey Pereverzev, **Toni Schneider**, Patrik Rorsman and Erik Renström (2005) Cav2.3 calcium channels control second phase insulin release. *J. Clin. Invest.* 115: 146-154.

Marco Weiergräber, Margit Henry, Michael Südkamp, Ernst-Rainer de Vivie, Jürgen Hescheler, **Toni Schneider** (2005) Ablation of Cav2.3 / E-type voltage-gated calcium channel results in cardiac arrhythmia and altered autonomic control within the murine cardiovascular system. *Basic Res. Cardiol.* 100: 1-13.

2004:

Klößner U., Pereverzev A, Leroy J, Krieger A, Vajna R, Pfitzer G, Hescheler J, Malécot CO, **Schneider T** (2004) The cytosolic II-III loop of Cav2.3 provides an essential determinant for the phorbol ester-mediated stimulation of E-type Ca₂₊ channel activity. *Eur. J. Neurosci.*, 19, 2659-2668.

Jahn Holger & **Schneider Toni** (2004) Diabetes mellitus als Komplikation der Behandlung mit atypischen Neuroleptika. Mögliche Pathomechanismen und Therapieempfehlungen. *Der Nervenarzt*, 75, 442-50.

Zhong-Ju Lu, Alexey Pereverzev, Hui-Lin Liu, Marco Weiergräber, Margit Henry, Andreas Krieger, Neil Smyth, Jürgen Hescheler, **Toni Schneider** (2004) Arrhythmia in Isolated Prenatal Hearts after Ablation of the Cav2.3 (a1E) Subunit of Voltage-gated Ca₂₊ Channels. *Cellul. Physiol. Biochem.*, 14, 11-22.

2003:

Jérôme Leroy, Alexey Pereverzev, Rolf Vajna, Ning Qin, Gabriele Pfitzer, Jürgen Hescheler, Claire Malécot, **Toni Schneider**, Udo Klößner (2003) A novel Ca₂₊-dependent regulation of neuronal E-type Ca₂₊ channels is mediated by an arginine-rich region in the cytosolic II-III linker. *Eur. J. Neurosci.*, 18, 841-855.

Atticus H. Hainsworth, Nicolle C McNaughton, Alexey Pereverzev, **Toni Schneider**, and Andrew D Randall (2003) Actions of sipatrigine, 202W92 and lamotrigine on R-type and T-type Ca₂₊ channel currents. *Eur. J. Pharmacol.* 467, 77 – 80.

Dirk Dietrich, Timo Kirschstein, Maria Kukley, Alexej Pereverzev, C. von der Brélie, **Toni Schneider**, Heinz Beck (2003) Functional specialization of presynaptic Cav2.3 Ca²⁺ channels. *Neuron*, 39, 483-496

2002:

Pereverzev Alexey, Jérôme Leroy, Claire Malécot, Gabriele Pfitzer, Jürgen Hescheler, Udo Klößner, **Toni Schneider** (2002) Alternate Splicing in the Cytosolic II-III Loop and Carboxy Terminus of the human E-type voltage-gated Ca₂₊ channels. - Electrophysiological characterization of isoforms. *Molec. Cellul. Neurosci.* 21:352-65.

Sochivko Dmitry, Alexey Pereverzev, Neil Smyth, Cornelia Gissel, **Toni Schneider**, Heinz Beck (2002) The a1E Ca₂₊ channel subunit underlies R-type Ca₂₊ currents in hippocampal and cortical pyramidal neurons. *J. Physiol.*, 542, 699-710.

Pereverzev Alexey, Marina Mikhna, Rolf Vajna, Cornelia Gissel, Margit Henry, Marco Weiergräber, Jürgen Hescheler, Neil Smyth, **Toni Schneider** (2002) Disturbances in glucose-tolerance, insulin-release and stress-

induced hyperglycemia upon disruption of the Cav2.3 (α1E) subunit of voltage-gated Ca²⁺ channels. *Mol. Endocrinol.*, 16, 884-95.

Pereverzev Alexey , Rolf Vajna, Gabriele Pfitzer, Jürgen Hescheler, Udo Klöckner, **Toni Schneider** (2002) Reduction of insulin secretion in the insulinoma cell line INS-1 by overexpression of a Ca²⁺ channel antisense-α1E cassette. *Eur. J. Endocrinol.*, 146, 881-889.

Michels Guido, Jan Matthes, Renate Handrock, Ute Kuchinke, Ferdi Groner, Leanne L. Cribbs, Alexey Pereverzev, **Toni Schneider**, Edward Perez-Reyes and Stefan Herzig (2002) Single-channel pharmacology of mibefradil in human native T-type and recombinant Cav3.2 calcium channels. *Mol. Pharmacol.*, 61, 682-694.

2001:

Vajna, R., Klöckner, U., Pereverzev, A., Weiergräber, M., Chen, X.-H., Miljanich, G., Klugbauer, N., Hescheler, J., Perez-Reyes, E., **Schneider, T.** (2001) Functional Coupling between 'R-Type' Ca²⁺ channels and insulin secretion in the insulinoma cell line INS-1. *Eur. J. Biochem.*, 268, 1066-1075.

2000:

Weiergräber Marco, Alexey Pereverzev, Rolf Vajna, Margit Henry, Martin Schramm, Wolfgang Nastainczyk, Heike Grabsch, **Toni Schneider** (2000) Immunodetection of α1E Voltage-gated Ca²⁺ Channel in Chromogranin-positive Muscle Cells of Rat Heart, and in Distal Tubules of Human Kidney. *J. Histochem. Cytochem.* 48, 807-819.

1999:

Klöckner, U., Lee, J.-H., Cribbs, L.L., Daud, A., Hescheler, J., Pereverzev, A., Perez-Reyes, E., **Schneider, T.** (1999) Comparison of the Ca²⁺ currents induced by expression of three cloned α1 subunits, α1G, α1H, and α1I of low-voltage-activated T-type Ca²⁺ channels. *Eur. J. Neurosci.* 11, 4171-4178.

Nakashima, Y.M., Pereverzev, A., **Schneider, T.**, Covey, D.F., Lingle, C.J. (1999) Blockade of Ba²⁺ current through human α1E channels by two steroid analogs, (+)-ACN and (+)-ECN. *Neuropharmacology* 38, 843-855.

Grabsch, H., Pereverzev, A., Weiergräber, M., Schramm, M., Henry, M., Vajna, R., Beattie, R., Volsen, S.G., Klöckner, U., Hescheler, J., **Schneider T.** (1999) Immunohistochemical Detection of α1E Voltage-gated Ca²⁺ Channel Isoforms in Cerebellum, INS-1 cells, and Neuroendocrine Cells of the Digestive System. *J.Histochem.Cytochem.* 47, 981-994.

Schramm, M., Vajna, R., Pereverzev, A., Tottene A., Tottene, A., Klöckner, U., Pietrobon, D., Hescheler, J., **Schneider, T.** (1999) Isoforms of α1E voltage-gated calcium channels in rat cerebellar granule cells. Detection of major calcium channel α1-transcripts by reverse transcription-polymerase chain reaction. *Neuroscience* 92, 565-575.

Jung-Ha Lee, Asif N. Daud, Leanne L. Cribbs, Antonio E. Lacerda, Alexei Pereverzev, Udo Klöckner, **Toni Schneider**, and Edward Perez-Reyes (1999) Cloning and expression of a novel member of the low-voltage-activated T-type calcium channel family. *J. Neurosci.*, 19, 1912-1921.

1998:

Nakashima, Y.M., Todorovic, S.M., Pereverzev, A., Hescheler, J., **Schneider, T.**, Lingle, C.J. (1998) Properties of Ba²⁺ currents arising from human α1E and α1Eβ3 constructs expressed in HEK293 cells: physiology, pharmacology, and comparison to native T-type Ba²⁺ currents. *Neuropharmacology*, 37, 957-972.

Vajna, R., Schramm, M., Pereverzev, A., Arnhold, S., Grabsch, H., Klöckner, U., Perez-Reyes E., Hescheler, J., **Schneider, T.** (1998) New Isoform of the Neuronal Ca²⁺ Channel α1E Subunit in Islets of Langerhans and Kidney. Distribution of Voltage-Gated Ca²⁺ Channel α1 Subunits in Cell Lines and Tissues. *Eur. J. Biochem.* 257, 274-285.

Pereverzev, A., Köckner, U., Henry, M., Grabsch, H., Vajna, R., Olyschläger, S., Viatchenko-Karpinski, S., Schröder, R., Hescheler, J., **Schneider T.** (1998) Structural Diversity of Voltage-Dependent Ca²⁺ Channel α1E-Subunit. *Eur. J. Neurosci.*, 10, 916-925.

1997:

Parent, L., **Schneider, T.**, Moore, C. P., Talwar, D. (1997) Subunit regulation of the human brain $\alpha 1E$ calcium channel. *J. Membrane Biology*, 160, 127 - 140.

Mehrke, G., Pereverzev, A., Grabsch, H., Hescheler, J., **Schneider T.** (1997) Receptor Mediated Modulation of Recombinant Neuronal Class E Calcium Channels. *FEBS Lett.* 408, 261-270

Schneider, T., Igelmund, P., Hescheler J. (1997) G Protein Interaction with K^+ and Ca^{2+} Channels. *Trends in Pharmacological Science*, 18, 8-11.

1995:

Perez-Reyes, E., **Schneider T.** (1995) The Molecular Biology of Calcium Channels. *Kidney International* 48, 111-1124.

Wei, X.Y., S. Pan, H. Kim, **T. Schneider**, E. Perez-Reyes & L. Birnbaumer (1995) Multiple subunits are required for the high affinity and allosteric regulation of dihydropyridine binding to the cardiac L-type Ca^{2+} channels. *J. Biol. Chem.* 270, 27106-27111.

Schneider T., Perez-Reyes, E., Nyormoi, O., Wei, X., Crawford, G.D., Smith R.G., Appel, S.H., Birnbaumer, L. (1995). Alpha-1 Subunits of Voltage Gated Ca^{2+} Channels in the Mesen-cephalon x Neuroblastoma Hybrid Cell Line MES23.5. *Neuroscience* 68, 479-485.

Timchenko, L., Nastainczyk, W., **Schneider, T.**, Hofmann, F., Caskey, C.T. (1995) Full-length myotonin protein kinase (72 kDa) displays serine kinase activity. *Proc. Natl. Acad. Sci. USA* 92, 5366 - 5370

1994:

Olcese, R., Qin, N., **Schneider, T.**, Neely, A., Wei, X., Birnbaumer, L., Stefani, E. (1994) The Amino Terminus of a Calcium Channel β Subunit Sets Rates of channel Inactivation Independently of the Subunit's Effect on Activation. *Neuron* 13, 1433-1438

Perez-Reyes, E., **Schneider T.** (1994) Ca^{2+} channels: Structure, function, and classification. *Drug Development Research* 33, 295-318

Nyormoi, O., **Schneider, T.**, Smith, R.G. (1994) A large scale preparative gel electrophoresis separation of $\alpha 1$ and $\alpha 2$ subunits of the voltage-gated Ca^{2+} channel from rabbit skeletal muscle. *Electrophoresis*, 15, 1183-1190.

Schneider T., X. Wei, Olcese, R., Costantin, J.L., Neely, A., Palade, P., Perez-Reyes, E., Qin, N., Zhou, J., Crawford, G.D., Smith, R.G., Appel, S.H., Stefani, E., Birnbaumer, L. (1994). Molecular analysis and functional expression of the human type E neuronal Ca^{2+} -channel $\alpha 1$ subunit. *Receptors and Channels* 2, 255 - 270.

Kimura, F., R. G. Smith, O. Delbono, O. Nyormoi, **T. Schneider**, W. Nastainczyk, F. Hofmann, E. Stefani & S. H. Appel (1994) Amyotrophic lateral sclerosis patient antibodies label the Ca^{2+} channel $\alpha 1$ subunit. *Ann. Neurol.* 35, 164-171.

1992 – 1989:

Smith, R.G., S. Hamilton, F. Hofmann, **T. Schneider**, W. Nastainczyk, L. Birnbaumer, E. Stefani & S. H. Appel (1992) Serum antibodies to L-type calcium channels in patients with amyotrophic lateral sclerosis. *N.Engl.J.Med.* 327, 1721-1728.

Regulla, S., **T. Schneider**, W. Nastainczyk, H. E. Meyer & F. Hofmann (1991) Identification of the site of interaction of the dihydropyridine channel blockers nitrendipine and azidopine with the calcium- channel $\alpha 1$ subunit. *EMBO J.* 10, 45-49.

Schneider, T., S. Regulla & F. Hofmann (1991) The devapamil-binding site of the purified skeletal muscle receptor for organic-calcium channel blockers is modulated by micromolar and millimolar concentrations of Ca^{2+} . *Eur.J.Biochem.* 200, 245-253.

Hofmann, F., V. Flockerzi, W. Nastainczyk, P. Ruth & **T. Schneider** (1990) The molecular structure and regulation of muscular calcium channels. *Curr.Top.Cell Regul.* 31, 223-239.

Hofmann, F., **T. Schneider**, A. Röhrkasten, W. Nastainczyk, M. Sieber, P. Ruth & V. Flockerzi (1989) Calcium channels: structure and function of receptors for calcium channel blockers in skeletal muscle. *Arzneimittelforschung.* 39, 164-168.

Hofmann F. & **Schneider T.** (1989) Calcium Channels: Structure and Function of the Skeletal Muscle Calcium Antagonist Receptor. In: Kotyk A. et al. (eds.) Highlights Modern Biochem. 1, 361-368

1988 – 1984: Jahn, H., W. Nastainczyk, A. Röhrkasten, **T. Schneider** & F. Hofmann (1988) Site-specific phosphorylation of the purified receptor for calcium-channel blockers by cAMP- and cGMP-dependent protein kinases, protein kinase C, calmodulin-dependent protein kinase II and casein kinase II. Eur.J.Biochem. 178, 535-542.

Hofmann, F., H. J. Oeken, **T. Schneider** & M. Sieber (1988) The biochemical properties of L-type calcium channels. J.Cardiovasc.Pharmacol. 12 Suppl 1, S25-S30.

Schneider, T. & F. Hofmann (1988) The bovine cardiac receptor for calcium channel blockers is a 195-kDa protein. Eur.J.Biochem. 174, 369-375.

Hofmann F., Nastainczyk W., Röhrkasten A., **Schneider T.**, Sieber M. (1987) Regulation of the L-type calcium channel. Trends in Pharmacol. Sci. 8, 393-398

Pönsngen-Schmidt E, **Schneider T.**, Hammer U, Betz A (1988) Comparison of Phosphoenolpyruvate-Carboxykinase from Autotrophically and Heterotrophically Grown Euglena and Its Role During Dark Anaerobiosis. Plant Physiology 86: 457-462.

Schneider T., Betz A (1985) Waxmonoester fermentation in Euglena gracilis T. Factors favouring the synthesis of odd-numbered fatty acids and alcohols. Planta 166: 67-73.

Schneider T., Borkowski C, Betz A (1984) Rotenone sensitive waxmonoester synthesis in anaerobic Euglena Gracilis. Evidence for operation of methylmalonyl-CoA pathway. In: Structure, Function and Metabolism of Plant Lipids (Siegenthaler P-A, Eichenberger W, eds), pp 577-580. United Kingdom: Elsevier Science Ltd.

Schneider T., Borkowski C, Betz A (1984) Wax Ester Formation in Euglena Gracilis During Anaerobiosis and Photoheterotrophic Growth. In: Advances in Photosynthesis Research. Proceedings of the VIth International Congress on Photosynthesis, August 1-6, 1983 (Sybesma C, ed), pp III.6.445-III.6.448. The Hague, The Netherlands: Martinus Nijhoff / Dr. W. Junk Publishers.